QuoreOne®

Holter e Looper cardíaco

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SAC
+ 55 11 973196913

quore.tech@quore.tech
www.quore.tech

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Every Heartbeat Matters.

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Cardiac Monitor and Loop Recorder

Instruction Manual

Contents

SAFETY INFORMATION	6
SYMBOLS USED	6
PROHIBITED ACTIONS	7
WARNINGS AND SAFETY	8
NOTES FOR ATTENTION	9
MESSAGES	11
CARE TO BE TAKEN	11
COMPATIBILITY AND ELECTRICAL SAFETY	11
INTRODUCTION	12
DESCRIPTION OF THE PRODUCT	12
PACKING WHEN SHIPPED	12
HOW DOES THE QUOREONE WORK?	12
INTENDED USE	13
CONTRAINDICATIONS	13
ACCOMPANYING MATERIAL	13
ACCESSORIES	13
TECHNICAL SPECIFICATIONS:	15
STATEMENT OF ESSENTIAL PERFORMANCE	15
STATEMENT ON APPLIED PARTS	17
CLASSIFICATION	17
PERFORMANCE CHARACTERISTICS	17
ELECTRICAL CHARACTERISTICS	18
PHYSICAL CHARACTERISTICS	18
ENVIRONMENTAL CHARACTERISTICS	18
GUIDELINES AND STATEMENTS ON ELECTROMAGNETIC EMISSIONS	19
INITIAL INSPECTION	22
HOW TO USE THE QUOREONE: INTRODUCTION	23
EXPECTED POSITION FOR USE OF THE QUOREONE	23
REPLACEMENT OF THE BATTERY	23
HAZARDOUS MATERIALS OR SUBSTANCES	24
Pause DEFINITION	24
CALCULATION OF HEART RATE	24
STARTING PROCEDURE	24
THE MULTIFUNCTION BUTTON FOR INTERFACE WITH PATIENT AND OPERATOR	25
MARKING AN EVENT DURING AN EXAMINATION	25
ALARMS AND INFORMATION SIGNALS	26
ADDITIONAL INFORMATION	27
DOWNLOADING AND INSTALLING THE QUOREAPP ON A SMARTPHONE	29

INSTALLING QUOREONE ON THE PATIENT	32
CONFIGURING THE QUOREONE	35
FUNCTIONS FREQUENTLY USED	42
REMOVING THE QUOREONE	42
SENDING THE DATA FROM QUOREONE TO THE QUOREAI PLATFORM.	43
GENERATION OF REPORT	44
CONFIGURATIONS OF THE QUOREONE EXAMINATION	44
<i>Duration of the examination</i>	45
<i>Channels</i>	45
RECHARGING THE BATTERY OF QUOREONE	46
SPECIAL PRECAUTIONS AND CARE TO BE TAKEN	47
SERVICE INFORMATION	47
<i>Requirements</i>	47
<i>Serial number</i>	47
<i>Cleaning, disinfecting and sterilization</i>	47
INFORMATION FOR SERVICE PERSONNEL	48
FIRMWARE VERSION	49
TERMINOLOGY AND DEFINITIONS	49
FCC HOMOLOGATION	50
CONTACT QUORETECH	50

CONTENTS – ILLUSTRATIONS

SAFETY INFORMATION	6
SYMBOLS USED	6
PROHIBITED ACTIONS	7
WARNINGS AND SAFETY	8
NOTES FOR ATTENTION	9
MESSAGES	11
CARE TO BE TAKEN	11
COMPATIBILITY AND ELECTRICAL SAFETY	11
INTRODUCTION	12
DESCRIPTION OF THE PRODUCT	12
PACKING WHEN SHIPPED	12
HOW DOES THE QUOREONE WORK?	12
INTENDED USE	13
CONTRAINDICATIONS	13
ACCOMPANYING MATERIAL	13
ACCESSORIES	13
TECHNICAL SPECIFICATIONS:	15
STATEMENT OF ESSENTIAL PERFORMANCE	15
STATEMENT ON APPLIED PARTS	17
CLASSIFICATION	17
PERFORMANCE CHARACTERISTICS	17
ELECTRICAL CHARACTERISTICS	18
PHYSICAL CHARACTERISTICS	18
ENVIRONMENTAL CHARACTERISTICS	18
GUIDELINES AND STATEMENTS ON ELECTROMAGNETIC EMISSIONS	19
INITIAL INSPECTION	22
HOW TO USE THE QUOREONE: INTRODUCTION	23
EXPECTED POSITION FOR USE OF THE QUOREONE	23
REPLACEMENT OF THE BATTERY	23
HAZARDOUS MATERIALS OR SUBSTANCES	24
PAUSE: DEFINITION	24
CALCULATION OF HEART RATE	24
STARTING PROCEDURE	24
THE MULTIFUNCTION BUTTON FOR INTERFACE WITH PATIENT AND OPERATOR	25
MARKING AN EVENT DURING AN EXAMINATION	25
ALARMS AND INFORMATION SIGNALS	26
ADDITIONAL INFORMATION	27
DOWNLOADING AND INSTALLING THE QUOREAPP ON A SMARTPHONE	29
INSTALLING QUOREONE ON THE PATIENT	32
CONFIGURING THE QUOREONE	35
FUNCTIONS FREQUENTLY USED	42
REMOVING THE QUOREONE	42
SENDING THE DATA FROM QUOREONE TO THE QUOREAI PLATFORM.	43
GENERATION OF REPORT	44
CONFIGURATIONS OF THE QUOREONE EXAMINATION	44

<i>Duration of the examination</i>	45
<i>Channels</i>	45
RECHARGING THE BATTERY OF QUOREONE	46
SPECIAL PRECAUTIONS AND CARE TO BE TAKEN	47
SERVICE INFORMATION	47
<i>Requirements</i>	47
<i>Serial number</i>	47
<i>Cleaning, disinfecting and sterilization</i>	47
INFORMATION FOR SERVICE PERSONNEL	48
FIRMWARE VERSION	49
TERMINOLOGY AND DEFINITIONS	49
FCC HOMOLOGATION	50
CONTACT QUORETECH	50

Safety Information

Symbols used

	<i>Date of manufacture</i>			<i>Consult the Instruction Manual</i>
	<i>Serial number</i>			<i>Warning</i>
	<i>Humidity limits for storage and transport</i>			<i>General warning symbol</i>
	<i>Temperature limits for storage and transport</i>			<i>Maximum stacking</i>
	<i>Keep away from rain</i>			<i>General symbol for "Prohibited action"</i>
	<i>Keep out of sunlight</i>			<i>General symbol for "Obligatory action"</i>
	<i>Fragile</i>			<i>Do not dispose of this device in common garbage</i>
	<i>Useful life</i>			<i>Applied Part: BF type</i>
	<i>This side up</i>		IP 44	<i>IP protection</i>

Prohibited actions

Never use QuoreOne during magnetic resonance procedures.

Never subject a patient wearing QuoreOne to an Automatic External Defibrillator (AED). Remove the QuoreOne before starting the procedure.

Do not use the QuoreOne on a patient with weight lower than 10 kg.

Never use the QuoreOne in any situation of travel or movement that has a three-point safety belt – this exposes the Patient to risk of serious injury, or even death.

Never submerge the QuoreOne in any liquid. Comply with the IP protection rating of QuoreOne.

If there is a need to submit the patient to the Heimlich Maneuver (removing obstruction to the upper airways), remove the QuoreOne before starting the procedure.

Cardiac Application: *This device cannot be used directly in contact with the heart.*

Never dispose of a QuoreOne or its batteries in common garbage. (Note: The internal battery is not accessible to the patient or the operator.)

NEVER perforate the QuoreOne to expose it to excessive heat or fire.
The QuoreOne has an internal lithium-ion battery that could explode, causing serious injury to health or even death of the Patient, Operator, Service personnel or any other person, or animal.

Operator: Never install a damaged QuoreOne on a patient. For all procedures of replacement or repair, contact Quoretech.

Warnings and safety

The equipment is not appropriate for use in an Oxygen-rich environment.

The equipment is not appropriate for use in the presence of flammable gases.

No modification or adaptation of this equipment is permitted.

Do not dismantle the equipment.

Patients with irritated or flaking skin should not use the QuoreOne.

Accidental spillage: To avoid electric shock or malfunction of the device, no liquid should be allowed to enter the device. If any liquid gets into the device, take the device out of use immediately and send it to technical assistance.

Leaks of electric current: Any failure of the equipment or failure in installation could result in electric shock to the patient.

Electrosurgery: If an electrosurgical device is used, the QuoreOne must be removed from the Patient.

Conducting materials: If an electrode comes into contact with conductive materials, this could cause electric shock or malfunction of the device. Keep the conducting parts of the electrodes and their related parts far from any other conducting items, including the earth (or an earth conductor).

Radio frequency interference: Known sources of radio frequency such as cell phones, radio or TV stations or other radio communication devices may cause adverse or unexpected operation of this device. For more information, consult specialized personnel.

Empty battery: If the battery runs out, immediately remove the QuoreOne device from the Patient and contact the Service personnel who installed it, reporting the event.

Notes for attention

It is advisable to obey the information on electromagnetic compatibility of this equipment described in this manual.

	<i>This product has no parts that can be adjusted or repaired by the Patient, an Operator or Service personnel. In the event of any doubt whatsoever, contact Quoretech.</i>
	<i>If there is an indication of low battery during use, the Operator should carry out the recharging procedure – see the section “Recharging the battery of QuoreOne”.</i>
	<i>When the “Electrode disconnected” Information Signal shows, follow the instructions described in the item “Installing the QuoreOne”.</i>
	<i>Safe use of this product is guaranteed to the standard of IEC 60601–1.</i>
	<i>The decision that this product is in operational condition is the responsibility of the Operator and their team. The Operator is responsible for determining that the equipment is safe to use, by carefully carrying out a pre-operational check. See the item (“Initial Inspection”).</i>
	<i>Read the whole of this manual before operating the equipment, and keep it for future reference and information.</i>
	<i>This product is not a pacemaker, nor a defibrillator, and does not have any functions similar to either of those devices.</i>
	<i>Functioning in the presence of pacemaker pulses: The QuoreOne is capable of registering the ECG signal in the presence of pulses from a pacemaker, but it is not capable of registering the activity of the pacemaker.</i>
	<i>The Operator is responsible for using the equipment in accordance with this Manual.</i>
	<i>This Manual is not intended as a substitute for appropriate and competent training. Operator: Apply to Quoretech for your training.</i>
	<i>After removing the equipment from its packaging, inspect it for any physical damage. If the equipment is damaged, do not install it on a Patient, but immediately contact Quoretech, giving the product’s serial number.</i>
	<i>Disposal material (including electrodes) is not supplied with the QuoreOne. This material should be acquired directly in the market. In the event of any doubts, contact Quoretech.</i>
	<i>An Internet connection is necessary for the QuoreApp to function. Internet service providers may charge for sending the data generated. Quoretech does not reimburse any such charges.</i>

Messages

There are safety messages at several points throughout this manual.

Care to be taken

Before starting use: Verify that the device is operating appropriately. See the chapter on “Service information”. If using the device with other devices or instruments, seek assistance from a specialist.

Disposal: This device is not disposable. At the end of its use, it should be returned to the same qualified health service that installed it. If you have any doubts about return of the device to the medical service, contact Quoretech or one of its authorized representatives.

Configuration of the equipment: The equipment must not be installed above, below, in contact with or close to any other equipment. If this situation cannot be avoided, it is necessary to carry out tests to check that the device is operating in the expected manner.

Interference: Electrical signals emitted by an electric blanket or thermal blanket may degrade the quality of the signal. Do not use this device jointly with an electric blanket or thermal blanket.

Instructions for use: For safe and continuous use of this equipment, the instructions listed must be followed. Note, however, that the instructions listed in this manual in no way replace established medical practices for attending a patient.

Compatibility and electrical safety

Take care: The QuoreOne requires special precautions in relation to EMC, and should be used in accordance with the information on EMC supplied in the item ***Electromagnetic immunity*** under ***Guidelines and Declarations*** in this manual.

Take care: Mobile radio frequency communication equipment may affect electromedical equipment.

Warning: The QuoreOne should not be used when stacked, nor when touching or resting on any other equipment.

Warning: The QuoreOne may suffer interference from other equipment, even if that equipment is in accordance with the requirements of CISPR.

Warning: Portable radio frequency communication equipment (including peripheral devices such as cables and antennae, internal or external) should be used at least 30 cm away from any part of the QuoreOne. If this distance is not complied with, the performance of QuoreOne may be degraded or even lost.

INTRODUCTION

Description of the product

QuoreOne is a Holter and a Loop Recorder. Its function is to acquire, store and transmit electrical signals generated by the sinus node of the cardiac muscle. Through the interpretation of the signals it is possible to diagnose disturbances in electrical generation and conduction and associate them with cardiopathies.

Packing when shipped

The product is shipped in a plastic cradle enclosed in a cardboard box enclosed in a plastic wrap.

Figure 1. The QuoreOne, packed.

How does the QuoreOne work?

Through electrodes (accessories) fixed on to the patient's skin, the QuoreOne acquires the electrical signals of the heart and stores them in an internal memory for subsequent wireless transfer of these signals to a digital platform – QuoreAI.

The technology:

QuoreOne uses low-energy Bluetooth wireless communication technology. Through the exclusive QuoreApp, a smartphone can be paired with the QuoreOne for the Operator to configure the device. After configuration, the QuoreOne will save all the signals acquired in non-volatile memory for subsequent transmission to the QuoreAI platform.

Note: The signals are transmitted to the platform exclusively by Quoretech service personnel: this process is not accessible to the Patient or the Operator.

Intended use

QuoreOne is designed for procedures known as Holter and Loop Monitoring. It acquires up to 3 channels of ECG electrical signals for up to 7 days – see the item “Performance characteristics”.

Contraindications

QuoreOne must not be used on any patient with weight less than 10 kg. For other restrictions on use, see the item “Prohibited actions”.

Accompanying material

The items shipped in the QuoreOne box are:

<i>Quantity</i>		<i>Code</i>		<i>Item</i>
<i>01</i>		<i>QO-01</i>		<i>QuoreOne</i>
<i>01</i>		<i>GR-01</i>		<i>QuoreOne Quick Guide</i>

Accessories

Accessory items – these are NOT shipped with the product:

The Instruction Manual is available only in digital version. To access and download the Manual visit the site of Quoretech. If you prefer you can enter the site on your smartphone by scanning the barcode below.

Figure 2. Download the QuoreOne Instruction Manual.

<i>Code</i>	<i>–</i>	<i>Item</i>
<i>EL-01</i>	<i>–</i>	<i>Monitoring electrode – foam (3M 2223 BRQ)</i>
<i>WP-01</i>	<i>–</i>	<i>Wireless induction charger</i>

Note: The accessory items EL-01 and WP-01 may be acquired directly in the market, without Quoretech as intermediary, provided the supplier is reliable and safe. To acquire them, use the specifications below.

Code	–	Description for acquisition in the market
<i>EL-01 foam monitoring electrode</i>	–	<i>Cardiac monitoring electrode for long-duration procedures, made of polyethylene foam, with hypoallergenic adhesive. Conducting adhesive gel, with metal pin in 304 or 316 stainless steel treated with silver/silver chloride.</i>
<i>WP-01</i>	–	<i>Wireless induction charger compatible with smartphones that use the Qi standard for wireless charging and Quick Charge.</i>

Recommended manufacturers and accessory models:

For EL-01	–	Item: Foam Monitoring Electrode Code No.: 2223BRQ Manufacturer: 3M
For WP-01	–	Item: Tabletop induction charger Code No.: QI15WA Manufacturer: Geonav

Use only wireless induction chargers certified by FCC.

Follow the information contained on the packaging of the electrodes. Never use an electrode that has been acquired with its packaging violated. Never use an electrode that is past its final use date or which has been inappropriately stored. Always obey the useful-life date of the electrode after the packaging has been opened. Non-compliance with these recommendations may result in loss of the product's Essential Performance.

In the event of any doubts, contact Quoretech.

TECHNICAL SPECIFICATIONS:

Statement of Essential Performance

The QuoreOne acquires, stores and transmits ECG signals for subsequent analysis of the data. When one or more electrodes has a bad contact or is disconnected, essential performance is lost and an Information Signal is emitted – see the item “Alarms and Information Signals”

Table 1. Requirements for Essential Performance

Requirement	–	Reference	–	Essential Performance Loss (Indications on the report / printed copy)
Heart rate	–	40 to 200 bpm	–	±5% deviation from measured heart rate and/or absence of / failure to highlight lower, average and higher heart rates, or the total number of heartbeats detected, on the printed report.
Supraventricular ectopy	–	Precocity is a factor established by the judgment of the doctor writing the medical report	–	±10% in relation to the factor set by the medical report author and/or absence of or failure to indicate totals for SVEB ¹ , or single SVEBs, paired SVEBs, measurements of TSV ¹ and duration of TSV (total beats), including the total number of each event, including an hour-by-hour table report, and the total.
Narrow ventricular ectopy	–	QRS with duration < 120ms	–	+12ms difference from Reference number, and/or absence or failure to indicate totals for VEB ¹ , single VEBs, paired VEBs, sequences of three or more VEBs or extension of runs (number of beats), and the totals of heartbeats for ventricular tachycardia of each episode. Absence of the number of minutes analyzed is also considered to be a Loss of Essential Function.
Wide ventricular ectopy	–	QRS with duration ≥ 120ms	–	12ms less than the Reference number
Bradycardia data	–	HR ≤ 50 bpm	–	Deviation of ± 3 bpm from Reference number and/or in the table report absence of, or failure to give, hour-by-hour totals of episodes of bradycardia, specifying the rate and duration of the episodes.

Table 2. Requirements for Essential Performance (continued)

Requirement	–	Reference	–	Essential Performance Loss – Evaluation on the report / printed copy
Pauses	–	RR interval ≥ 2.5 seconds.	–	± 200ms from Reference number and/or absence of indication or failure to indicate the total number of pauses detected, and location and duration of the longest pauses.
Shifts of the ST segment	–	Not available	–	Not applicable
Printed copy of the ECG	–	Print with well-defined charts and obligatory fields shown.	–	Failure in generation of the file for printing. Failure or absence from the ECG tape of 25mm/s and/or indication of channel. Failure to show, or absence of, the time on the tape, heartbeat frequency and/or annotation of the tape. Failure in or absence of identification of the Patient on the tapes (or at least one per page). Generation of file for printing in a format not A4 or Letter. Failure in or absence of channel calibration signal on each tape.

- 1 –
 SVEB – supraventricular ectopic heartbeat
 VEB – ventricular ectopic heartbeat
 SVT – supraventricular tachycardia
 HR – heart rate

Statement on Applied Parts

QuoreOne has contact with the patient through the electrodes (see the item “Accessories”). The other parts of QuoreOne are considered to be Accessible Parts.

Classification

Type of protection against electric shock	–	The equipment is energized internally.
Degree of protection against electric shock	–	This is low frequency equipment. The voltages of the internal power source are sufficiently low for there not to be a perception of electric shock under any conditions.
IP protection	–	IP44
Degree of protection for use in the presence of inflammable anesthetic mixed with air or with nitrous oxide	–	This equipment is not appropriate for use in the presence of inflammable anesthetics mixed with air or with oxygen or nitrous oxide.
Degree of protection for use in the presence of an Oxygen-rich environment.	–	This equipment is not appropriate for use in an Oxygen-rich environment.
Method of sterilization or disinfection recommended by the manufacturer.	–	This product should not be submitted to any sterilization process. For cleaning procedures consult the item “Cleaning, disinfection and sterilization”.
Modes of operation	–	Continuous operation
Is the equipment transportable?	–	The equipment is portable.

Performance characteristics

ECG channels	–	1, 2 or 3 channels. Selectable.
Memory capacity	–	Up to 7 days
Recording format	–	Continuous
Time of use	–	Up to 7 days
Storage time	–	6 months
Storage time outside packaging	–	Up to 30 days
Time to recharge internal battery	–	Approximately 2h30min (depending on health of the battery and model of charger.)

Electrical characteristics

Type of Applied Part	–	Low frequency (electrode – accessory)
ECG response frequency	–	0.67 to 40 Hz
ECG input impedance	–	$\geq 10 \text{ M}\Omega$
ECG differential range	–	$\pm 1.65 \text{ mV}$
ECG A/D sampling rate	–	200 Hz
ECG resolution	–	16 bits
RF transmission / reception	–	2.4 GHz Low energy Bluetooth Effective RF power < 1mW
Transmission band frequency	–	2.4 GHz

Physical characteristics

Dimensions	–	130 x 140 x 23 mm
Weight (with battery and without packaging)	–	60 grams

Environmental characteristics

Operating temperature	–	5°C to 43 °C ¹⁻²
Operating altitude	–	-1 to 4,000 meters
Humidity for operation and storage	–	10% to 95% (without condensation)
Temperature for transport and storage (short-term, rapid transport)	–	10°C to 35 °C
Temperature for long-term storage	–	10°C to 35 °C
Storage altitude	–	-1 to 4,000 meters
IP Classification	–	IP44
Expected useful life	–	5 years

Notes:

¹ During recharging of the battery, the QuoreOne may reach temperatures higher than operational temperature – read the item “Battery recharging procedure”

² During an examination, the QuoreOne may occasionally reach temperatures higher than 41°C, but will never go above 43°C. When a temperature above 41°C occurs, the patient will feel a slight discomfort due to the heat generated at the fixing points of the QuoreOne on the thorax. If temperature exceeds 43°C, the QuoreOne should be immediately removed from the Patient due to risk of burns and/or damage to the Patient’s tissue. If the sensation of heat is uncomfortable, it is probably because part of the main body of the QuoreOne is touching the patient’s skin. In this case, it is recommended that the Patient should reposition their posture and members so that only the electrodes have contact with the thorax, leaving space for ventilation between the body of the QuoreOne and the patient’s skin. Older patients usually have more sensitive skin and may have a

stronger sensation of discomfort. Avoid using any topical medication under or in an area near the QuoreOne. The reaction to such medications may be stronger due to the high temperature. In normal conditions of use, the QuoreOne does not exert any surface pressure on the Patient. If the posture of the Patient's body causes a superficial pressure to be exerted between the product and the patient the sensation of heat may be stronger. In these cases, reorient the posture to relieve the pressure exerted. In the event of any doubt about the heat generated and its consequences, contact Quoretech.

Guidelines and statements on electromagnetic emissions

This device is compliant with FCC regulations. Operation is subject to the two conditions:

- (1) This device may cause damaging interference; and
- (2) This device is susceptible to any interference received, including interference that may cause undesired operation. This device does not have accessories. The use of any accessory may possibly not ensure compliance with FCC's guidelines for exposure to radio frequency.

Any changes or modifications not expressly approved by the party responsible for compliance may annul the user's authority to operate the equipment. The QuoreOne has been tested and complies with FCC's RF exposure levels when used and operated for the purpose for which it is intended, and as instructed in the Manual.

Table 3. Guidelines and statements on electromagnetic emissions

QuoreOne is intended for use in the electromagnetic environment specified in the tables below. The user of QuoreOne must ensure that it is used in such an environment.		
Emission Test	Compliance	Orientation on electromagnetic environment
CISPR 11 RF emissions	Group 1	QuoreOne uses a radio frequency system only for communication via Bluetooth technology with a cell phone. As a result, the levels of RF emission are very low, and it is not likely that any interference will be caused to any electronic equipment close to it. The RF system of QuoreOne is active only during its configuration. The RF system of the QuoreOne is not active during execution of the examination.
CISPR 11 RF emissions	Class B	QuoreOne is appropriate for use in all establishments, including domestic establishments.
IEC 61000-3-2 harmonics emissions	Not applicable	Not applicable
Voltage fluctuation / IEC scintillation emission	Not applicable	Not applicable

Table 4. Guidelines and statements on electromagnetic immunity

<p>QuoreOne is intended for use in the electromagnetic environment specified in the table below. The client or user of QuoreOne must ensure it is used in an environment which has the restrictions listed below.</p>			
Immunity test	Test level	Compliance level	Orientation on the electromagnetic environment
Electrostatic discharge (ESD) – IEC 61000-4-2	± 8 kV by contact ± 15 kV via air	± 8 kV by contact ± 15 kV via air	The floor should be of wood, concrete, or ceramic. If the floor is covered with a synthetic material, relative humidity should be at least 30%.
Immunity test	IEC 60601 test level	IEC 60601 compliance level	Orientation on the electromagnetic environment
Electrostatic discharge (ESD) – IEC 61000-4-2	± 8 kV by contact ± 15 kV via air	± 8 kV by contact ± 15 kV via air	The floor should be of wood, concrete, or ceramic. If the floor is covered with a synthetic material, relative humidity should be at least 30%.
Immunity test	IEC 60601 test level	IEC 60601 compliance level	Orientation on the electromagnetic environment
Radiated RF: IEC 61000-4-3	<p>10 V/m 80 MHz to 2.7 GHz</p> <p>28 V/m 385, 450, 810, 870, 930 MHz 18 Hz pulsed 9 V/m 710, 745, 780 MHz 217 Hz Pulsed</p> <p>28 V/m 1720, 1845, 1970, 2450 MHz 217 Hz Pulsed</p> <p>9 V/m 5240, 5500, 5783 MHz 217 Hz Pulsed</p>	<p>10 V/m</p> <p>28 V/m</p> <p>9 V/m</p> <p>28 V/m</p> <p>9 V/m</p>	<p>Portable or mobile RF communication equipment, including cables, should not be used close to a QuoreOne. The recommendation for separation distance is calculated by this equation applied to the transmission frequency: Recommended separation distance: $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ 80 MHz at 800 MHz $d = 2.3\sqrt{P}$ 800 MHz at 2.5 GHz – where P is the maximum output power of the transmitter in Watts (W) as informed by the transmitter manufacturer, and <i>d</i> is the recommended separation in meters.</p> <p>The intensities of the fields of fixed RF transmitters, as determined by a test of the local electromagnetic field ^a, should be lower than the compliance level in each frequency range ^b. Interference may occur in equipment that is close to QuoreOne and has the symbol:</p> 

Network frequency (50/60 Hz) Magnetic field IEC 61000-4-8	30 A/m	30 A/m	Magnetic fields coming from the electricity network should be at levels characteristic of a typical commercial or hospital environment.
---	--------	--------	---

NOTE 1 – At 80 MHz to 800 MHz, the higher frequency interval is applied.

NOTE 2 – These guidelines may possibly not apply in all situations. Propagation of electromagnetic waves is affected by absorption by and reflection from structures, objects and people.

^a It is not possible to make any precise theoretical prediction of the magnetic fields from fixed stations such as radio base stations (cell phone/wireless), mobile radios, amateur radio, AM or FM radio stations or TV transmitters.
To evaluate the electromagnetic environment caused by the presence of RF transmitters, an electromagnetic test survey should be considered.
If the strength of the field measured at the location where the QuoreOne is used exceeds the RF compliance level, the QuoreOne should be observed to check that it is functioning normally.
If abnormal performance is detected, additional measures such as reorienting or repositioning the QuoreOne may be necessary.

^b Above the frequency band from 150 kHz to 80 MHz, the intensity of the field should be lower than 3 V/m.

Table 5. Recommended separation distance between the QuoreOne and any portable or mobile RF communication equipment.

The QuoreOne is designed for use in an electromagnetic environment where interference from RF disturbances is controlled. The client or user of QuoreOne can help to avoid electromagnetic interference by maintaining a minimum distance between the portable RF communication equipment (transmitter) and the QuoreOne, as recommended in the table below, in accordance with the maximum output power of the communication equipment.

Maximum output power of the transmitter (watts)	Separation distance for different transmitter frequency bands		
	150 kHz to 80 MHz	80 MHz to 800 MHz	80 MHz to 2.5 GHz
	$d = 1.2\sqrt{P}$	$d = 1.2\sqrt{P}$	$d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmission sources with maximum output power not listed above, the recommended separation distance in meters (m) can be determined using the equation applied to the transmitter frequency, where P is the maximum exit potential in watts (w) as stated by the transmitter manufacturer.

Note 1. In the range from 80 MHz to 800 MHz, the separation distance for the highest frequencies is applied.

Note 2. It is possible that these guidelines may not apply in all situations.

Electromagnetic wave propagation is affected by absorption by and reflection from structures, objects and people.


INITIAL INSPECTION

When receiving the QuoreOne, check that the packaging is unbroken and sealed. If there are any signs of violation of the packaging, refuse receipt of the QuoreOne or contact your supplier to arrange replacement. If any issues or questions remain, contact Quoretech.

Check for any sign of penetration of water or dust inside the box of the QuoreOne. Refuse to use it if it is wet or dirty.

The packaging must contain the 'Rapid Guide' with technical information and instructions on how to download the Instruction Manual. Note: This manual is not supplied in printed form together with the QuoreOne. You can make a printed version of this Manual if necessary.

Check whether the QuoreOne is complete and unharmed, and that its parts – including the Multifunction Button and connectors of the electrodes (Figures 3 and 4) are included.

The QuoreOne has an internal battery that is rechargeable only by induction technology. It does not have any connectors or cables for recharging. To acquire an induction charger, please see the item "Accessories" in this Manual.

Check that the serial number, and the other marks on the product, are complete. Refuse to use it if any of the marks are illegible.

HOW TO USE THE QUOREONE: INTRODUCTION

Expected position for use of the QuoreOne

The Patient may use the QuoreOne in the standing, sitting or lying down position.

Situations that should be avoided when the QuoreOne is installed:

- Avoid sports or any other activities that produce abundant sweat, because sweat may result in reduced adherence of the electrode to the skin.
- When entering any environment with a temperature considered to be extreme, such as a sauna, cold store, etc., obey the temperature and humidity limits for use stated in this Manual, and the specified degree of protection for the QuoreOne.

The expected position for the procedures of installation, configuration and removal, both for the Patient and for the Operator, are described in the respective items of this Manual.

For use of the QuoreOne in any situations other than those specified in this Manual, contact Quoretech before beginning to install the QuoreOne.

Replacement of the battery

The battery of the QuoreOne is internal, and rechargeable. It can be replaced only by Quoretech. When it is close to the end of its useful life, Quoretech will arrange exchange of battery.

If you identify any problems related to low battery autonomy, contact Quoretech.

Patient, Operator and Service personnel: Do not try to dismantle the QuoreOne.

Service personnel: There is no fuse or component requiring adjustment in the interior of the QuoreOne.

The battery of QuoreOne is recharged exclusively by an induction charger. See the item "Accessories". The QuoreOne does not have any connectors or cables for recharging the battery.

Hazardous materials or substances

If there is a leak from the battery, remove the device immediately, and immediately inform Quoretech. If the leak touches the skin, remove the electrodes and wash the skin with abundant soap and water.

Ideally, the device and electrode should be removed only by a qualified and trained professional, but if there is a leak of the battery, it is understood that the risk of the patient removing the device is less than that of the patient remaining with the QuoreOne when it is in failure mode. For this, follow the steps described in the item "Removing the

QuoreOne”, requesting a second person to remove the device from the Patient. The Patient should not remove the device himself/herself, except when no other option is available.

Pause: definition

For the purposes of operation of QuoreOne, an R-R interval of 2.5 seconds or more is defined as a Pause.

Calculation of heart rate

Heart rate is calculated by this equation:

$$HR = 1,000 \times \frac{60}{R-R}$$

Where:

HR = heartbeat rate in beats per minute;

RR = the RR interval of the prior heartbeat cycle in milliseconds.

Starting procedure

To turn on the QuoreOne, press the Multifunction Button (Figure 3) for seven seconds. The Status Indicator light will turn on, indicating that the QuoreOne has been successfully started.

QuoreOne has only one control button, with the functions of (i) turning on, (ii) turning off, and (iii) marking an Event. Consult the item “Multifunction button for interface with Patient and Operator” for more details on its working.

For fuller interpretation of the Status Indicator light, see the item “Status Light and Information Signals”.

Read the item “Installing the QuoreOne on the Patient” for the information needed for the Operator to start operation of the QuoreOne, including configuration using the QuoreApp.

The Multifunction Button for interface with Patient and Operator

The QuoreOne has only one button, with multiple function, providing the interface with the patient and the operator. It carries out different functions depending on the status of the QuoreOne at the time.

Status of QuoreOne	–	Press the Multifunction Button for this period of time	–	Function executed
Turned off	–	7 seconds	–	Turns on the QuoreOne

Turned on, but not carrying out examination	–	7 seconds	–	Turns off the QuoreOne
Turned on, and examination in progress	–	Any period of time	–	Marks an Event

It is not possible to turn off the QuoreOne while an examination is in progress.

Marking an Event during an examination

Correlation between the time of occurrence of an external event (for example, administration of drugs, presence of symptoms, physical activities) and the ECG complexes of patients is essential for clinical interpretation.

QuoreOne has an Event Marking function. If during an examination the patient wishes to insert marking of an event, pressing the multifunction button for one second will make this mark visible to the Operator, or the doctor preparing an opinion, when analyzing the signals captured by the examination.

Up to 10 events can be marked per hour of examination, with minimum intervals of 10 seconds between each marking.

Pressing the Multifunction button does not create a request for help or support. If necessary, ask any person close by for help, or ring the emergency service.

Alarms and Information Signals

The QuoreOne does not have an alarm.

The QuoreOne has two information signals: (i) low battery, and (ii) electrode disconnected or with bad contact.

Low Battery:

When QuoreOne detects a low level of charge in the battery, it will inform the Patient or Operator through the Status Light (see “Status Light and Information Signals”).

If the battery is less than 90% full, QuoreOne will not allow an examination to start.

Electrode disconnected, or with bad contact:

When QuoreOne detects one or more electrodes with a bad contact, or disconnected, it will inform the patient or operator in two ways: Through the Status Light, and also through the QuoreApp – the latter only during configurations of an examination.

Status Light and Information Signals

The situations below establish a relation between the Status Light and the Information Signals corresponding to each one of the situations of the QuoreOne.

When the QuoreOne is turned on using the Multifunction Button:

- The **green light** remains on during the phase of initialization, and turns off when that phase is finalized.
- The **green light** flashes **twice** after the Multifunction Button has been pressed for seven seconds, and QuoreOne is turned off.

When QuoreOne is turned on, and is not carrying out an examination:

- The **green light** flashes **every five seconds**, indicating state of readiness.
- The **green light** flashes **three times** when Bluetooth is paired.
- The **green light** flashes **five times** when the examination is about to begin.
- The **blue light** flashes **twice** every five seconds to indicate that the battery is low.

When QuoreOne is re-charging the battery on the induction charger:

- The **green light** flashes **three times** when QuoreOne recognizes the wireless induction charger.
- The **blue light** flashes **twice every five seconds**, indicating that QuoreOne is being recharged.
- The **blue light** remains on **continuously**, indicating that QuoreOne is completely charged.

When QuoreOne is carrying out an examination:

- The **green light flashes** when the multifunction button is pressed, indicating that an Event has been marked.
- The **green light** flashes **ten times**, indicating that the examination has been finalized and that QuoreOne will turn itself off.
- The **blue light** flashes **once every two seconds**, indicating that one or more electrodes are disconnected or have bad contact.
- The **blue light** flashes **three times every five seconds**, indicating that the battery is low and that one or more of the electrodes is disconnected or has bad contact.
- The **blue light** flashes **three times every five seconds**, indicating that the battery is low and that one or more of the electrodes is disconnected or has bad contact.

Note. When carrying out an examination and the battery level is good, and all electrodes are properly connected, the Status Light remains off. If the Patient or Operator has any doubt about the situation of the QuoreOne device (e.g. whether it is turned on, and/or carrying out an examination), press the Multifunction Button for one second and see whether the Status Light lights up **green**. If it does, it indicates that QuoreOne is turned on and that an examination is in progress. It also **lights up** to inform the user that an Event has been marked. Take note of the time when this check was made, and inform the operator for this Event Mark to be left out of the Report.

Additional information

QuoreOne is fed exclusively by battery, and has no connection with any outside source of electricity. There is no voltage equalization point or conductor.

QuoreOne does not have connection cables. Electrodes (see the item “Accessories”) must be connected directly to the connection points of QuoreOne. Although the arms of the QuoreOne are flexible, neither the terminals nor the electrodes should ever be placed in contact with each other.

Degree of protection: QuoreOne has IP44 protection:

- protection against solid objects with a diameter of 1 mm or more; and
- Protection against water splashing from all directions.

QuoreOne, therefore, should not be used in an environment with excessive dust, due to the risk of penetration of the product, causing interference in the reading or even irreversible damage to the device.

Patient: Never use QuoreOne during activities such as: entering a swimming pool, an Ofuro bath, a river, the sea, a sauna, or any other similar environment.

QuoreOne does not have a function to pause an examination. QuoreOne detects when one or more channels are disconnected from the patient, such as, for example, when electrodes are changed.

If it is necessary to remove QuoreOne from the patient for any activity, such as, for example, travel by car, it is not necessary to take off the electrodes, nor any action to pause the examination. The time during which QuoreOne was without signal will not create an impact for the final results. The algorithms of the artificial intelligence platform will remove this period of absence of signals from the final report.

QuoreOne was not designed to detect and/or measure movements in the ST segment, thus analyses of this segment are not carried out.

Downloading and installing the QuoreApp on a smartphone

QuoreApp is a smartphone app used to configure QuoreOne.

The Operator has to have access to a smart phone for QuoreApp to be installed.

These instructions are directed to the Operator who will install require the QuoreApp application on the smartphone.

The QuoreApp should be downloaded ONLY from Google Play Store. Do not download QuoreApp from any other source.

Installing QuoreApp requires a smartphone with Android operating system (Version 8.0 or later).

*QuoreApp is **NOT** compatible with iOS (Apple iPhone), Microsoft or Blackberry.*

*QuoreApp is **NOT** compatible with tablets, notebooks, personal computers or similar devices.*

An Internet connection is necessary to download QuoreApp. Data transmitted may be charged by your Internet service provider. Quoretech does not reimburse any such charges.

Download the QuoreApp™ application from Google Play Store. For this, you will need to search for QuoreApp in Google Play Store, or use the camera of your Smartphone to read the QR code – See Figure 5.

Figure 5. Link for downloading the QuoreApp.

After installing QuoreApp, open it on your smartphone.

It will request permission to access your phone (Figure 6) and your files (Figure 7). In both situations, click on the option “ALLOW”. If you have pressed the “DENY” option one or more times, uninstall the app and reinstall it.


Figure 6. QuoreApp requesting permission to access a smartphone.


Figure 7. QuoreApp requesting access to the files of a smartphone.

After the permissions have been given, the screen in Figure 8 will open. The information on this screen refers to the smartphone.

- **“Wi-Fi connection”**: Indicates when your telephone is duly paired with QuoreOne.
- **“Battery”**: Indicates the level of charge of the battery of the smartphone.
- **“Bluetooth”**: Indicates when your telephone is duly paired with QuoreOne.
- **“Storage:”** Indicates how much memory of your smartphone is available to store, among other data, the signals collected by QuoreOne.
- **“Exams in queue for sending”**: Tells you whether there are examinations stored and ready to be sent.

- **“Server Connection”**: Informs the status of the QuoreApp’s connection with the data server.


Figure 8. QuoreApp installed, and showing information on the smartphone.

OK! Your QuoreApp is now duly installed.

Installing QuoreOne on the Patient

These instructions are for the Operator who will install QuoreOne on the patient. If you are a Patient, do not try to install QuoreOne on yourself. Seek a professional to install it on you.

QuoreOne should only be installed by people trained and qualified for these activities.

Installation of QuoreOne by an untrained person may cause damage to the Patient's skin, or damage the product.

Never use QuoreOne on patients whose weight is less than 10kg.

If the Patient has an injury, red skin, or itching on the region to be shaved, suspend installation of QuoreOne until the skin is healthy. Ask the patient whether the region to be shaved and cleaned is sensitive.

Material necessary for installation of QuoreOne:

<i>Item</i>	<i>Quantity</i>
QuoreOne device	1
Monitoring electrodes (see the item "Accessories")	4
Razor	1
Disposable gauze compress	1
Isopropyl alcohol solution	25 ml

Expected position of the Operator for installation of QuoreOne:

The professional who will install QuoreOne on the Patient should position her/himself in front of the patient at a distance of between 30 and 40 cm.

Procedure for shaving, cleaning, fixing of the electrodes and installation of QuoreOne:

Note: Keep in mind that the instructions below are fundamental for ensuring good quality in the signal acquired.

Use the razor to shave the region indicated in Figure 9. Considering the line connecting the two nipples (red line in diagram), shave the area indicated by the rectangles. Note that two-thirds of the height of the vertical rectangle is above the line connecting the nipples. Note: Although some patients do not have apparent hair on this region of their body, it is recommended that the procedure should be carried out in any event, regardless of the patient's race, gender, weight, height or age.



Figure 9. Area to be shaved.

With a disposable gauze compress soaked in isopropyl alcohol solution, clean the shaved area, by pressing the compress vigorously against the patient's skin. This procedure should last around three minutes. If necessary, replace the gauze during the procedure, always using gauze soaked with isopropyl alcohol solution.

Separate one of the foam electrodes (Figure 10) and remove the plastic protection film (Figure 11).

Install the electrode on the patient's thorax in the position illustrated in Figure 12. With your fingertips, vigorously press the electrode against the thorax of the patient to ensure good adhesion. Install the other three electrodes (without removing the plastic protection film) on the QuoreOne, as shown in Figures 13 and 14.

Bring the QuoreOne close to the thorax, holding it by its main body – Figure 15.

Fit the central connector of the QuoreOne (highlighted in Figure 14) on the metallic terminal of the electrode already installed on the patient. Press the QuoreOne against the patient's thorax until a click is heard, indicating that the connection has been successful.

Figure 12. Installation of the central electrode.

Figure 13. Fitting the electrode to QuoreOne.

Figure 14. QuoreOne with the electrodes fitted. The arrow indicates the central electrode.

Figure 15. Bringing the QuoreOne close to the thorax for installation.

Figure 16. QuoreOne with the central electrode connected to the electrode on the body.

After QuoreOne has been fitted to the central electrode and is suspended on the patient's thorax (Figure 16), remove the protective film from the electrode closest to the left arm (Figure 17), and, aligning QuoreOne, pressure it, by the extremity only, against the patient's thorax – Figure 18. Repeat this with the electrode closest to the right arm, and finally with the lower electrode.

If any one of the electrodes does not achieve the appropriate adhesion, remove it and repeat the process of shaving, cleaning and installation of a new electrode.

Figure 17. Removing the protective film from the electrode.

Figure 18. Pressuring the electrode to ensure good quality of adhesion.

Figure 19. The QuoreOne duly installed.

Configuring the QuoreOne

These instructions are addressed to the Operator, who will configure the QuoreOne using the QuoreApp.

The Operator must be the only person with access to configuration of QuoreOne. The Patient does not have access to the QuoreApp application.

The following messages refer exclusively to the smartphone that is in use in the process of configuring QuoreOne.

To configure the QuoreOne, the Operator must have the QuoreApp application installed. See the item “Downloading and installing the QuoreApp on the smartphone” in this Manual.

Before starting the configuration of QuoreOne, close all the other applications that are being used.

Do not share the Bluetooth connection with any other process than QuoreOne.

Do not attend calls, whether from the cellular mobile network or through apps, during configuration of QuoreOne.

Expected position for configuration of QuoreOne:

The Operator should remain within a radius of 1m from the patient for configurations via Bluetooth to be effective.

Configuration procedure:

Open QuoreOne on your smartphone. The opening screen, Figure 20, will indicate the current status of the smartphone.

The information and warnings below refer to the opening screen.

An active Internet connection is necessary for configuration of QuoreOne.

QuoreOne will not start execution of an examination if the battery level of the smartphone is less than 90%.

Do NOT start configuration of QuoreOne if the battery level of the smartphone is lower than 20%.

The Bluetooth of your smartphone must be active and not paired with any other device.

*If your smartphone is connected to the Internet and the status of the connection with the server is not “OK”, do not start any configurations or examination. Immediately contact Quoretech technical support. **This message applies to the smartphone, both of the Operator and of the Patient.***

On the opening page, insert your User name and Token.

The ‘Token’ is a password supplied by a separate application which generates Tokens. Operator: Before starting your first access to the QuoreApp, contact Quoretech (contact details are at the end of this Manual) for instructions for client registry and authorization of access to the token generator.


Figure 20. Data of the smartphone.


Figure 21. Input the User name and Token.

On the screen showing registered patients – Figure 22 – select the Patient who will be submitted to examination with QuoreOne. If the patient is not yet listed, click on the option “NEW EXAM”.


Figure 22. Patients registered.

On the screen for registering a new examination, fill in the patient’s details and click “SAVE” – Figures 23 and 24.


Figure 23. Insert (“New examination”).

Figure 24. Additional information for new examination.

On the “Synchronize device” screen, Figure 25, look for the QuoreOne to be installed on the list of devices available for pairing. Click on the “+” corresponding to the QuoreOne to be paired, for the pairing process to start.

The QuoreOne must be turned on for it to be found on QuoreApp. To turn the QuoreOne on, consult the item “Starting Procedure”.


Figure 25. QuoreOne devices available for pairing.


Figure 26. QuoreOne paired.


Figure 27. Configuration of examination.


Figure 28. Indication of installation of a QuoreOne.

When the pairing has been successfully completed, the indication “Paired” will appear on the QuoreApp screen – Figure 26.

On the “Configurations” screen, select the duration of the examination (in minutes, hours or days) and how many channels should be monitored – Figure 27. Also see the item “Functions frequently used”.

After the configuration is completed, click “Next”. QuoreApp will indicate that the QuoreOne is ready to test signal quality. Confirm whether the QuoreOne is duly installed on the patient

- see the item “Installing QuoreOne on the patient”, and click on the option “Test QuoreOne”
- Figure 28.

The test of QuoreOne will be started – Figures 29 and 30.


Figure 29. starting test...


Figure 30. Collecting test data.

Wait until the screen with ECG signals is shown – Figure 31.

Observe the quality of the signals acquired and the connection status of the electrodes of the Holter – see the item “Alarms and Information Signals”. If you want to repeat the test of the Holter, click on the option “go back” on the taskbar of your smartphone and repeat the test of the Holter, as shown in Figure 29.

Click on “Next” to finalize configuration, and wait for the request for confirmation of the Test Signal Quality to appear.

If the test signal quality is considered adequate (This assessment is for the Operator’s judgment), press “YES”. If not, press “NO” to return to the test of the signal strength. If the signal quality is not yet adequate, check the electrodes, and if necessary replace them.

Figure 31. Visualization of signal and status of the electrodes.


Figure 32. Confirmation of signal quality.


Figure 33. Configuration completed.


Figure 34. Loading list of patients.

When the message "Configuration completed" appears, this means that configuration has ended and the examination has started successfully – Figure 33. After checking the information of the test that has been started, press the option "Confirm". At this moment, the QuoreApp prepares to exhibit the list of the names of patients already registered for execution of a new configuration procedure – Figure 34.

The Patient now has the QuoreOne duly installed and may be released for his/her daily activities, subject to the restrictions stated in this Manual.

Functions frequently used

Holter procedures are frequently used for a period of 24 hours with three channels.

Removing the QuoreOne

These instructions are for the Operator, who will remove the QuoreOne from the patient. If you are a Patient, do not try to remove the QuoreOne yourself. Get a professional to remove it for you.

When to remove the QuoreOne?

The QuoreOne should be removed from the patient after completion of the period of the examination.

Note 1. It is not necessary to turn QuoreOne off, because it is automatically turned off when the period of examination ends.

Removal of the QuoreOne from the Patient must be carried out only by people trained and qualified for these activities.

Removal of the QuoreOne by an untrained person may cause injury to the Patient's skin, and/or damage the product.

Materials necessary for removal of QuoreOne from the patient:

<i>Item</i>	<i>Quantity</i>
Disposable gauze compress	1
isopropyl alcohol solution	15ml

Expected position for removal of the QuoreOne:

The professional who will remove the QuoreOne from the patient should position him/herself in front of the patient at a distance of between 30 and 40 cm.

Procedure for removal of the electrodes and cleaning of the QuoreOne:

Disconnect the extremities of the QuoreOne from the electrodes, pulling by the plastic parts of the QuoreOne as indicated in Figure 35.

Pull the QuoreOne by its main body, holding it between your thumb and first finger (Figure 36). After removing the QuoreOne, the electrodes will still be installed on the Patient's thorax (Figure 37).


Figure 35. Removing the QuoreOne device from the electrode.

Figure 36. Removing the QuoreOne from the patient's thorax.

With fingertips, remove all the electrodes from the patient's thorax, pulling only lightly so that there is no damage to tissue (Figure 38).

Throw the electrodes away, in the normal garbage disposal.

With a disposable gauze compress soaked in isopropyl alcohol solution, clean the area of the electrode, pressing the compress softly against the patient's skin. Repeat this until there is no residue of the glue from the electrodes on the patient's thorax.

If the Patient's skin shows signs of irritation, orient the patient to seek medical assistance for analysis and treatment.


Figure 37. The Patient's thorax after removal of the QuoreOne device.

Figure 38. Remove the electrodes slowly with fingertips.

Sending the data from QuoreOne to the QuoreAI platform.

Operator: Uploading of the examination data is an activity exclusively carried out by the technical team of Quoretech. The QuoreApp does not have a function for data loading. For transfer of the data, send the QuoreOne to Quoretech for the data to be loaded and inserted in the QuoreAI platform.

QuoreApp does not have a function of uploading of the examination carried out.

The Patient, the Operator and Service personnel do not have access to the data uploading procedure.

Generation of report

Generation of the final report or medical Opinion is the responsibility of Quoretech.

Finalizing the examination, Quoretech will generate the report – through the QuoreAI platform – and this will be made available through the QuoreApp to the Service personnel who configured the QuoreOne on the patient.

Quoretech has a team of doctors to write reports, who fully analyze the signals captured, classify them, and produce the report based on the signals generated.

The report is always generated using the web browser Mozilla Firefox (version evaluated: 74.0). Tests made with the browser Google Chrome show low resolution of the report, especially of the gridlines, and this can adversely affect evaluation of the ECG signal.

Configurations of the QuoreOne examination

Configuration of QuoreOne comprises two items:

- (i) Duration of the Examination; and
- (ii) Number of Channels – Figure 39.

For instructions for pairing Bluetooth and other means of communication, consult the item “Procedures for configuration” and Figure 25.


Figure 39. Configuration screen.

Duration of the examination

QuoreOne can carry out examinations with a duration from one minute to up to 7 days uninterrupted (see the intervals available on the selector of the QuoreApp).

On the configuration screen, the horizontal bar has a scale in number of minutes, hours or days, for choosing the desired duration of the examination. The bar and the green sphere indicate the Duration chosen for the Examination.

Channels

The QuoreOne device has three different ‘derivations’ or channels of monitoring. The central electrode, located in the central part of the device, is known as the reference electrode, and it is not possible to choose between its being activated or not. Select which channel you want to be activated (by pressing, for example, “Channel 1”), for the examination to be started. The channels selected for monitoring will be presented in **green**; while those not monitored will be shown in **black**.

For example: If all the channels are **green**, it indicates that the examination is being carried out with three channels.

Differences between the channels: The three channels of the QuoreOne producing different measurement characteristics, with the following polarities:

C1	–	Terminals RA (+) and LA (–)
C2	–	Terminals RA (+) and LL (–)
C3	–	Terminals LL (+) and LA (–)

Where: RA – Right arm; LA – Left arm; LL – Left leg.

The channels used are shown on each ECG tape in the examination report.

Recharging the battery of QuoreOne

This procedure applies only to the Operator or Service Personnel:

Patient: Never, when an examination is in progress, place the QuoreOne on top of any wireless induction charger – due to the risk of the examination being interrupted or stopping, and loss of essential performance of the QuoreOne.

QuoreOne has an internal rechargeable lithium-ion battery. Recharging the battery requires a wireless induction charger, which is not supplied. See the item “Accessories” in this Manual.

Use only wireless induction chargers certified by FCC.

To recharge the battery: position the wireless induction charger on a non-metallic, flat, clean, dry table with no metal object within a radius of 20 cm. Check the compatibility of the working voltage of the charger and that of the nearest plug and connect it. Avoid use of extensions, voltage transformers, or multiple connectors of any kind.

Read the instructions in the user manual of the wireless induction charger.

Place the QuoreOne delicately on the wireless induction charger, with the Status Light and Multifunction Button facing upwards. Position the QuoreOne on the charger so that the Multifunction Button aligns on the vertical axis with the geometrical center of the wireless induction charger, or if the positioning is indicated on the charger, position the Multifunction Button in the center of the area shown – Figure 40.

As soon as the QuoreOne is positioned on the charger, the Status Light will indicate that battery charging has started. Consult information on the Status Light in the section “Status Light and Information Signals” in this manual

QuoreOne does not have a visual indicator of the level of battery charge.

The QuoreApp does not indicate the level of battery charge during the recharging process.

Figure 40. QuoreOne being recharged.

SPECIAL PRECAUTIONS AND CARE TO BE TAKEN

Service information

Requirements

Consult the Quoretech maintenance service through the contacts shown in the item: “CONTACT Quoretech”. Any attempt at repair by any non-authorized technical service will cancel the equipment’s guarantee. It is the responsibility of the user to report any need for service to Quoretech or to one of its authorized agents.

Serial number

The device’s serial number is unique, and can be used to identify it. The serial number is visible on the back part of the device.

Cleaning, disinfecting and sterilization

NEVER subject the QuoreOne to any procedure of sterilization.

After each use, the QuoreOne should be cleaned and disinfected with a soft flannel slightly humidified with 70% INPM hydrated ethyl alcohol solution.

If you see or suspect a battery leak, send the QuoreOne for maintenance. Do not dismantle the QuoreOne.

The procedure of cleaning and disinfection should be carried out in a clean, well-aired, well-illuminated area, with good air circulation at ambient temperature and at ambient temperature. There is no set limit to the number of times the QuoreOne can be submitted to the procedure of cleaning and disinfection, but the legibility of its markings should be checked with each new cleaning. For more information on the effect of multiple cleaning and the state of conservation of the QuoreOne, consult the item. “Information for Service Personnel” in this Instruction Manual.

Use latex surgical gloves, a mask that covers at least the nose and mouth, and protection goggles to carry out this procedure.

Information for Service Personnel

Minimum qualifications: A “Service Person” is defined as a professional with a minimum training at technical level in electronics or telecommunications, with training specified under ‘general items’ of IEC 60601-1, ISO 14971, and with adequate and competent training carried out and registered by Quoretech on “General Aspects of the QuoreOne”.

- The QuoreOne does not have any internal fuse.
- The QuoreOne’s internal power supply is not replaceable.
- There is no internal component that can be adjusted.
- To replace the battery, send the QuoreOne to Quoretech.

Calibration: QuoreOne should be calibrated periodically, the period being defined as one year from the first use since the last calibration carried out. For calibration procedures, see the item “Statement” below.

Statement: Quoretech does not make available any circuit diagrams, list of components, descriptions or instructions for calibration. Any calibration, adjustment or repair, when necessary, should be carried out exclusively by Quoretech. For this, contact us through the channels indicated at the end of this instruction manual.

Check the following – at least once every 10 examinations:

- i. Legibility of the markings of the QuoreOne. Execution of multiple cleaning and disinfection may affect the legibility of these markings. Observe and monitor wear over the useful life of the QuoreOne, and when necessary send the QuoreOne to Quoretech before the markings become illegible.¹
- ii. Good condition of the parts of the QuoreOne.
 - a) The connectors of the metal electrodes should be free of oxidation and with no signs of crashing, or any other damage.
 - b) Rigid parts must be complete, without cracks or fissures, peeling or wear.
 - c) The flexible parts (the arms) must have their flexibilities preserved, without drying out, alteration of color or any cuts or tears.

If one or more of the above items is not compliant, send the QuoreOne to Quoretech for maintenance.

- 1- Since the upload of data is carried out by the Quoretech team, a finely detailed assessment of the condition of conservation and wear of the product will be carried out with every 100 examinations realized. The monitoring of this interval is carried out by the Quoretech team. Whenever necessary, parts may be replaced to maintain the completeness of the product.

Firmware version

Firmware 1.0.18

Terminology and definitions

Statement: The following terminology is taken from IEC 60601-1 2010 and its 2016 update. Some information may have been suppressed only for the purpose of summarizing the text, but this shall be without prejudice to interpretation in its application to QuoreOne.

Patient: A living person or animal subject to a medical, surgical or dental procedure. A Patient can be an Operator.

Operator: A person who works with the equipment. *In the case of QuoreOne, a Patient is exclusively the person on whom the QuoreOne is installed, and this person does not exercise the function of operator.*

Operator: A person who works with the equipment. *In the case of QuoreOne, the operator is the professional who installs and configures the QuoreOne on the patient.*

Service Personnel: The individuals or entity that is legally or morally responsible for the Responsible Organization that installs, assembles, maintains and repairs Electromedical Equipment.

Procedure: A specific way of executing an activity or process.

Process: Group of interrelated or interactive activities that transforms inputs into products or outputs.

Accessible part: Any part of an item of electrical equipment, different from an Applied Part, which may be touched with a **standard test finger**.

Applied Part: Any part of an item of electromedical equipment which in normal use comes into physical contact with the patient for the electromedical equipment to carry out its function.

Basic safety: Absence of Unacceptable Risk directly caused by physical dangers when the electromedical equipment is used in normal conditions, or in abnormal condition under only one single failure.

Continuous operation: Operation in normal use for an unlimited time without the specified temperature limits being exceeded.

Equipment envelope: Exterior surface of the electrical equipment or its parts.

Essential performance: Performance of a clinical function not related to Basic Safety, in which loss or degradation of the limits specified by the manufacturer results in an Unacceptable Risk.

Expected useful life: Period of time specified by the manufacturer during which it is expected that the Electromedical Equipment or Electromedical System remains safe for use (that is to say, maintains Basic Safety and Essential Performance).

Internally energized: This term refers to electrical equipment that is able to operate from an internal electrical power supply.

Internal electrical power supply: A source of electrical power for operation of the equipment, which is part of the equipment and produces an electric current from another form of energy. Examples: Chemical, mechanical, solar or nuclear.

Normal condition: Condition in which all the available means for protection from dangers remain intact.

Intended Use, or Intended Purpose: Use of the product, Process or service in accordance with the specifications, instructions and information supplied by the manufacturer.

Normal use: Operation, including inspection routine, regulation by any Operator, and state of readiness, in accordance with the instructions for use.

Oxygen-rich environment: An environment in which the concentration is greater than 25% under environmental pressure up to 110 kPa, or, for environmental pressure exceeding 110 kPa, where the partial pressure of oxygen is greater than 27.5 kPa .

Accessory: Additional part for use with the equipment so as to: provide the Intended Use; adapt it to some special use; facilitate its use; improve its performance; or enable its functions to be integrated with those of other equipment.

Other terminology used in this Manual:

QuoreApp: An **application** for smartphone developed by Quoretech to configure and receive the data from the QuoreOne.

QuoreAI: A digital platform of Quoretech which receives the data transmitted through the QuoreApp, and processes and generates reports.

CONTACT QUORETECH

Quoretech LLC

Technical responsible: Eng. Marcos Vinicius Pinto e Silva - CREA-SP 5061757301

200 Bellevue Parkway Suite 210
Wilmington
DE 19809
USA
+1 206 6270030

Every Heartbeat Matters.

quore.tech@quore.tech

www.quore.tech

Designed by Quoretech in Brazil

Made in Brazil

Version of this manual : Rev. 4 EN